



AMERIGLIDE UNIVERSAL 400 PR



AmeriGlide's Universal 400 PR adds Powered Rotation (PR) to the Universal 400 framework, which makes it even easier to lift your power chair or scooter in and out of your automobile.

With its 400 pound lifting capacity, this power chair and scooter lift can hold nearly any power chair/scooter on the market. Besides being a powerful lift, the Universal 400 PR features: an integrated offset post that ensures that your power chair or scooter lifts clear of your vehicle's bumper, a base that adjusts to better match your vehicle's dimensions, and a manually adjustable telescoping boom (with manually adjustable length) that makes it very easy to attach the lifting strap to your electric wheelchair or adult scooter.

The Universal 400 PR has a very pleasing look that will match most auto interiors and comes standard with a remote hand control, weather resistant exterior, and a removable hook.

Safely stow your power chair or scooter inside your vehicle and out of nature's way!

Specifications

| | |
|--------------------------------------|--------------------------|
| Lifting Capacity | 400lbs |
| Operation | Power Lifting & Rotation |
| Adj Boom Height | 30" – 40" |
| Adj Boom Length | 21" – 32" |
| Pitch | 5 degrees |
| Adj Rotation Limit | Yes |
| Weather Resistant | Yes |
| Docking Device Included | Yes |
| Passenger Side Standard | Yes |
| Driver Side Upgrade Available | Yes |
| Warranty: | 3-Years |

