



Bariatric Rollator with 8-inch Wheels

The **ProBasics Bariatric Aluminum Rollator** is made with an extra-strength, reinforced aluminum frame and has large reinforced forks on 8-inch wheels for improved maneuverability and stability. The padded seat and removable padded backrest provide a comfortable resting area. It features easy-to-operate ergonomic hand brakes, a water-resistant nylon pouch, and a built-in fold strap for easy storage and transport. Supports patient weights up to 400 pounds



Extra Wide Seat and Water-Resistant Pouch



Easy-to-Operate Ergonomic Hand Brakes



Reinforced forks and 8-inch Wheels

FEATURES

- Lightweight aluminum frame
- Two-button removable curved padded backrest
- Comfortable, extra wide flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 8" black casters with reinforced forks and soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0149
E0156 (seat)

COLORS

Burgundy
Finish



Blue
Finish



SPECIFICATIONS

Overall Width:	29"
Width b/w Hand Grips:	23.5"
Seat Dimensions:	19" (W) x 13" (D)
Seat-to-Floor Height:	22"
Weight Capacity:	400 lb
Adjustable Handle Height:	33.5" – 38.5"
Wheel Size:	8"
Product Weight:	18.7 lb
Carton Dimensions:	29.1" x 8.6" x 36.2"
Shipping Weight:	22.7 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLAB8BG	ProBasics Bariatric Aluminum Rollator	Burgundy	815067071746	1/cs
RLAB8BL	ProBasics Bariatric Aluminum Rollator	Blue	815067071739	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.