



## Buzzaround GP130



The new Buzzaround CarryOn folding power wheelchair is truly in a class all by itself! This unique power chair is super portable with a 21" turning radius and center wheel drive for optimal coverage and maneuverability making it the perfect partner for your next adventure. Our exclusive PowerLink™ Connectivity makes assembly and disassembly a snap! Independent 4-wheel suspension and a springless articulating chassis provide a smooth ride indoors and out. The CarryOn features a 300 lb. weight capacity, front and rear light package, and textured soft touch armrests. Enjoy long range coverage up to 12.5 miles on the airline friendly 12.5 AH lithium-ion battery, which can be removed and taken as a carry on while traveling by air, sea or train. The heaviest piece is 46 lbs. (without seat and batteries) and the height when folded is only 10" for small space storage tucked out of the way. Carry on with life in the CarryOn power chair!



### Standard Colors:



Velvet Red  
(Installed)



Coastal Blue  
(Included in Box)

### Optional Colors:



Crisp White



Canary Yellow




Stone Grey



Orchid Pink



Army Green  
(Matte)



Desert Tan

**Only 46 Lbs.!\***

\*Without Seat & Battery

**Most Portable, Foldable,  
Lightweight Power Chair with  
Center Wheel Drive!**

**Tightest Turning Radius!**

### Standard Features:

- Center Wheel Drive for Optimal Control & Maneuverability
- Independent 4-Wheel Suspension
- Springless Articulating Chassis
- USB Port
- Front & Rear Lights
- Width Adjustable Armrests
- Textured Soft Touch Armrests
- Angle Adjustable Foot Plate
- Easy Access, Airline-Friendly, 12.5 AH Lithium-ion Battery with Travel Case
- PowerLink™ Connectivity for Easy Disassembly/Assembly
- Removable Storage Basket
- Easily Interchangeable Color Shrouds

# BUZZ AROUND CarryOn

Power Wheelchair

## FOLDING POWER WHEELCHAIR

<b>Model Number</b>	<b>GP130</b>
Medicare Code	N/A
Weight Capacity	300 lbs.
Drive Wheels	Center
Maximum Speed	3.85 mph
Operating Range <sup>1</sup>	12.5 miles
Ground Clearance	2"
Turning Radius	21"
Type Batteries	12.5 AH Lithium-Ion
Colors	Red, Blue
Freewheel Mode	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board, 2A
Controller	40 Amp
<b>Dimensions:</b>	
Length	38"
Width	24"
Height (ground to top of back of seat)	34.25"
Ground to Top of Seat	21"
<b>Weight of Unit: (assembled)</b>	<b>74 lbs.</b>
Base	46 lbs.
Seat	22 lbs
Battery	6 lbs.
<b>Tires:</b>	
Drive Wheels	7.75"
Casters	5.5"
<b>Seat:</b>	
Back Height - no headrest	15.5"
Width x Depth	18" x 17"
Color	Black
<b>Warranty:</b>	
Frame	5 Years
Drive Train	13 Months
Electronics <sup>2</sup>	13 Months
Batteries	12 Months

<sup>1</sup> Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain, and battery type.  
<sup>2</sup> Electronics warranty excludes batteries.



Battery & battery travel bag standard



Folds to only 10" high for storing in small spaces

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Phone: 800-624-6374  
 Fax: 800-628-5165  
 www.goldentech.com

 Find us on Facebook  
[www.facebook.com/goldentechnologies](http://www.facebook.com/goldentechnologies)

SHT-GP130-08234